

imaet PRO Version 11.2 Release Notes

Software

New Features & Improvements

- Health History Form on Client Profile restored, functional, and savable.
- New Disclaimer Pop-Up Window on Log-In page
- New Injury & Trauma category of frequencies added to the Athletic Panel

Box and Accessories

- None

Minimum System Requirements

- Windows 10 or later
- 8 Gigabytes RAM
- i5 Processor
- 512 Gigabytes of Internal Memory

Known Issues & Limitations

- The Biofield Panel is intended for advanced users. Newer users should utilize the Allergen Panel instead until having completed QBA training or after a minimum of 6 months of use.
- Panel Info pages are unavailable on the Pet Panel and Frequency Specific Panel
- A glitch on the Allergen Panel can occur when a user does not click the Clear button to remove the data in the Search Field and instead overwrites the previous search term
- The Bio-Scan on the Geopathic Panel continues to display after scan is completed. Does not affect functionality
- Anti-virus or anti-malware software running during use can cause malfunction
- Pending Windows Updates and driver updates can cause malfunction
- Opening the program while it is already running, creating multiple instances of the program, causes malfunction.

See our Recommended [Device Best Practices](#) for optimal performance